

The design of our models is protected against counterfeiting.<sup>TS</sup>

## Seat belt

*Belt to prevent people from falling down or falling over in a chair or seat.  
These are made of a comfortable material and can be washed at 60°C and air-dried.  
Quality 1065 : Breathable, flame-retardant anti-tear material.*



### **Md 31920 D<sup>TS</sup>**

Seat belt with an extra Close-up® zipper to attach a lap blanket.

### **Md 31920 DC<sup>TS</sup>**

Lap blanket in microfleece, dimensions 80/85 cm, zipable.

**Md 31920 DB<sup>TS</sup>**: Lap blanket in qual. 1010 dim. 80/85 cm, zipable.

#### *Purpose:*

- prevent the patient from sinking into the wheelchair or care seat
- provide support so that the patient can adopt a good sitting posture
- increase patient safety and prevent the patient from falling by immobilising in the wheelchair or care seat.
- To be used for mentally confused persons, who are not aggressive.

#### *Quality:*

1065, permanently flame retardant, breathable material,  
Washable at 60°C and tumble-dryer proof. Short drying time at low temperature.

#### Washing instructions for quality 1010 and 1065:

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* DON'T ADD FABRIC CONDITIONER.
3. *Dry* Dry short at low temperature in the drier or dry hanging.  
Don't dry together with hydrophilic materials (such as cotton,...).