

Nursing aids

The design of our models is protected against counterfeiting. TS

Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.



Md 40631 TS

Cover for the boomerang pillow in quality 1010 with press studs over the underside.

Washing instructions for guality 1010:	
1. Washing	Wash separately at 60°C with <u>little</u> detergent.
2. Rinse	DO NOT ADD FABRIC CONDITIONER.
3. Dry	Dry <u>very briefly</u> at <u>low</u> temperature in the dryer or hang to dry.
	Do not dry together with moisture-retaining textiles (like cotton,).
1	

