

The design of our models is protected against counterfeiting.<sup>TS</sup>

## Seat belt

*Belt to prevent people from falling down or falling over in a chair or seat.  
These are made of a comfortable material and can be washed at 60°C and air-dried.  
Quality 1065 : Breathable, flame-retardant anti-tear material.*



### Md 31920<sup>TS</sup>

Seat belt blue, prevents sliding down, one size, fastening on the back of the seat by means of a buckle which must be opened with two hands.



#### Purpose:

- prevent the patient from sinking into the wheelchair or care seat
- provide support so that the patient can adopt a good sitting posture
- increase patient safety and prevent the patient from falling by immobilising in the wheelchair or care seat.
- To be used for mentally confused persons, who are not aggressive.

**Quality:** 1065, permanently flame retardant, breathable material,  
Washable at 60°C and tumble-dryer proof. Short drying time at low temperature.

#### Washing instructions for quality 1010 and 1065:

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* DON'T ADD FABRIC CONDITIONER.
3. *Dry* Dry short at low temperature in the drier or dry hanging.  
Don't dry together with hydrophilic materials (such as cotton,...).