

**Nursing aids** 

The design of our models is protected against counterfeiting. TS

## **Positioning cushions**

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.



## Md 40656 TS

Cover for the straight support cushion, length +/- 60 cm, in quality 1010 with zipper.

Washing instructions for quality 1010:1. WashingWash separately at 60°C with little detergent.2. RinseDO NOT ADD FABRIC CONDITIONER.3. DryDry short at low temperature in the drier or dry hanging.<br/>Don't dry together with hydrophilic materials (such as cotton,...).

