

Nursing aids

The design of our models is protected against counterfeiting. TS

Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.



Md 40656 TS

Cover for the straight support cushion, length +/- 60 cm, in quality 1010 with zipper.

Washing instructions for quality 1010:1. WashingWash separately at 60°C with little detergent.2. RinseDO NOT ADD FABRIC CONDITIONER.3. DryDry very briefly at low temperature in the dryer or hang to dry.
Do not dry together with moisture-retaining textiles (like cotton,...).

