

Nursing aids

The design of our models is protected against counterfeiting. TS

Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.



Md 40686 TS

Cover for the abduction cushion, in quality 1010 with zipper.

Washing instructions for quality 1010:1. WashingWash separately at 60°C with little detergent.2. RinseDO NOT ADD FABRIC CONDITIONER.3. DryDry short at low temperature in the drier or dry hanging.
Don't dry together with hydrophilic materials (such as cotton,...).

