

*The design of our models is protected against counterfeiting. TS*

## Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

*Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.*



### Md 40651<sup>TS</sup>

Cover for the support cushion for paralysis, outer circumference +/- 285 cm, in quality 1010 with zipper.

#### Washing instructions for quality 1010:

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* DO NOT ADD FABRIC CONDITIONER.
3. *Dry* Dry very briefly at low temperature in the dryer or hang to dry.  
Do not dry together with moisture-retaining textiles (like cotton,...).

