

The design of our models is protected against counterfeiting.<sup>TS</sup>

## Side rail protector

The side rail protector is used to prevent restless residents from injuring themselves at the side rail.

Quality 1010: Breathable, flame-retardant textile with high tear resistance.



### Md 40058 <sup>TS</sup>

Side rail protector with flame retardant foam, polyurethane cover. Surrounded by a cover in 1010 with buckles for fitting around the side rail. This model is only protecting around the side rail.

#### Purpose:

- Prevent the patient from hurting himself or protruding limbs through the siderail.
- To be used for mentally confused persons who are not aggressive.

#### Quality:

\*1010, permanently fire retardant, breathable material, drying cabinet resistant.  
Short drying time at low temperature.  
\*PU, wipeable

#### Washing instructions for quality 1010 :

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* DON'T ADD FABRIC CONDITIONER.
3. *Dry* Dry short at low temperature in the drier or dry hanging. Don't dry together with hydrophilic materials (such as cotton,...).