

The design of our models is protected against counterfeiting.^{TS}

Side rail protector

The side rail protector is used to prevent restless residents from injuring themselves at the side rail.

Quality 1010: Breathable, moisture-repellent and flame-retardant textile with high tear resistance.



Md 40058^{TS}

Side rail protector with flame retardant foam and polyurethane cover. Surrounded by a cover in 1010 with buckles for fitting around the side rail. This model is only protecting around the side rail.

Purpose:

- Prevent the patient from hurting himself or protruding limbs through the siderail.
- To be used for mentally confused persons who are not aggressive.

Quality:

*1010, permanently flame-retardant, breathable material, drying cabinet resistant.
Short drying time at low temperature.
*PU, wipeable

Washing instructions for quality 1010 :

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* DON'T ADD FABRIC CONDITIONER.
3. *Dry* Dry very briefly at low temperature in the dryer or hang to dry. Don't dry together with moisture-retaining textiles (like cotton,...).