

The design of our models is protected against counterfeiting. ^{TS}

Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.



Md 40646 ^{TS}

Cover for the neck cushion, outer circumference +/- 100 cm, in quality 1010 with zipper.

Washing instructions for quality 1010:

1. Washing Wash separately at 60°C with little detergent.
2. Rinse DO NOT ADD FABRIC CONDITIONER.
3. Dry Dry very briefly at low temperature in the dryer or hang to dry.
Do not dry together with moisture-retaining textiles (like cotton,...).