

Anti-tear linen and restraint

The design of our models is protected against counterfeiting. TS

Anti-tear clothing and bed linen (1065)

For patients with severe disposition to tearing we offer a specific range of products, such as: anti tear quilts, sheets, shirts, short and T-shirt,.... Quality 1065: Breathable, moisture-repellent and flame-retardant anti-tear material.



Md 40201 TS

Anti-tear sheet 250/140 cm.

Washing instructions for quality 1065:

1. Washing Wash separately at 60°C with <u>little</u> detergent.
2. Rinse DO NOT ADD FABRIC CONDITIONER.

3. Dry Dry very briefly at low temperature in the dryer or hang to dry.

Do not dry together with moisture-retaining textiles (like cotton,...).

