

Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

Quality 1010: Breathable, moisture-resistant, and flame-retardant fabric with high tear resistance.



Md 40631 TS

Cover for the boomerang pillow, circumference +/- 200 cm, in quality 1010 with press studs over the underside.

Washing instructions for quality 1010:

1. *Washing* Wash separately at 60°C with little detergent.
2. *Rinse* **DO NOT ADD FABRIC CONDITIONER.**
3. *Dry* Dry very briefly at low temperature in the dryer or hang to dry.
Do not dry together with moisture-retaining textiles (like cotton,...).

