

## Positioning cushions

Positioning cushions are used to support parts of the body or the whole body and to place the body in the desired position.

*Quality 1010 : Breathable, moisture-repellent and flame-retardant fabric with high tear resistance.*



### Md 40681<sup>TS</sup>

Cover for the chair cushion with 3 compartments, in quality 1010 with zipper.

#### Washing instructions for quality 1010:

- 1. Washing* Wash separately at 60°C with little detergent.
- 2. Rinse* DO NOT ADD FABRIC CONDITIONER.
- 3. Dry* Dry very briefly at low temperature in the dryer or hang to dry.  
Do not dry together with moisture-retaining textiles (like cotton,...).

