

*The design of our models is protected against counterfeiting.<sup>TS</sup>*

## **Anti-tear clothing and bed linen (1065)**

*For patients with severe **disposition to tearing** we offer a specific range of products, such as: anti tear quilts, sheets, shirts, short and T-shirt,....*

*Quality 1065 : Breathable, moisture-repellent and flame-retardant anti-tear material.*



### **Md 40131<sup>TS</sup>**

Anti-tear pillow, washable 60/70 cm.

Washing instructions for quality 1065:

- |                   |   |
|-------------------|---|
| 1. <i>Washing</i> | Wash separately at 60°C with <u>little</u> detergent.   |
| 2. <i>Rinse</i>   | <u>DO NOT ADD FABRIC CONDITIONER.</u>   |
| 3. <i>Dry</i>     | Dry <u>very briefly</u> at <u>low</u> temperature in the dryer or hang to dry.<br>Do not dry together with moisture-retaining textiles (like cotton,...). |